

SIS^{plus} WHEN CONTROL COUNTS

When you're faced with deep tissue trauma, whether due to injury, disease or surgery, you want a product that will allow you to control the healing process as much as possible. SIS^{plus}, a unique SIS product processed with ion therapy, can be used as a surgical patch for tissue repair and reinforcement.


What It Is...

- Porcine-derived small intestine submucosa (SIS)
- Processed to maintain a framework of scaffolding for tissue ingrowth
- Contains structural proteins (such as collagen and elastin) as well as other extracellular matrix proteins and GAG compounds

What It Does...

- Supports and promotes wound healing and tissue remodeling by providing a biocompatible, absorbable tissue repair matrix
- Contains angiogenic promoters, including a unique combination of ions which promote rapid healing and a reduced infection potential
- Offers excellent physical properties, including high suture retention and bursting strength
- Encourages rapid tissue ingrowth
- Mitogenic
- Anti-microbial

SIS^{plus} offers you all the benefits you've come to expect from an SIS product, along with the increased rapid healing response that comes with our unique ion therapy.



Actual Size

Product	Price
Single Laminate (10cm x 7cm)	\$120
Multi Laminate (4 ply; 10cm x 7cm)	\$200
Ocular Disc (4/pkg.; 10mm)	\$60
Ocular Disc (4/pkg.; 15mm)	\$60
<ul style="list-style-type: none"> • Supplied sterile in a dry configuration • For custom sizing, please call 	